

Standard Operating Procedure

Reagent Handling, Labeling & Cold-Chain Storage Control

Document No.	Revision	Effective Date	Review Cycle
SOP-QC-014	4.0	15 June 2026	24 months
Department	Site	Owner	Status
Quality Assurance	Nexa Biosciences – Meridian Campus	L. Fontaine, QA Manager	✔ Approved for Use

This document is controlled. Printed copies are uncontrolled unless stamped in red ink by the Document Control Office. Verify this is the current revision at the Nexa QMS portal before use in any GxP activity.

1 Purpose and Scope

1.1 Purpose

This Standard Operating Procedure (SOP) defines the controlled process for receiving, labeling, storing, and monitoring laboratory reagents at Nexa Biosciences, including temperature-sensitive and light-sensitive materials. The objective is to prevent reagent degradation, cross-contamination, and traceability gaps that could invalidate downstream analytical or cell-based assay results.

1.2 Scope

This procedure applies to all reagents, buffers, standards, and consumable chemical stocks received into the Meridian Campus wet-lab suites (Rooms 2A, 2B, and 3C), including:

- ▶ Commercially sourced reagents requiring 2–8 °C, –20 °C, or –80 °C storage.
- ▶ In-house prepared buffers and working solutions assigned an internal lot number.
- ▶ Reference standards used for analytical method qualification and QC release testing.

This SOP does **not** cover controlled substances (see SOP-QC-021, *Controlled Substance Custody*) or biological samples under chain-of-custody (see SOP-QC-009).

1.3 Regulatory Basis

This procedure is written to satisfy the storage-and-handling expectations of ISO/IEC 17025:2017 §6.4 and the internal Nexa Quality Manual QM-01, and is cross-referenced during annual internal audits [1].

Deviations from any numbered step in Section 4 must be logged per SOP-QC-002 (*Deviation and CAPA Reporting*) before the affected reagent lot is released for use.

2 Roles and Responsibilities

Role	Responsibility
Laboratory Technician	Executes receiving, labeling, and storage steps; records lot data in the reagent log within 1 hour of receipt.
Laboratory Supervisor	Verifies weekly cold-chain logs; co-signs storage-excursion investigations; approves reagent disposal requests.
QA Reviewer	Audits reagent log entries monthly; confirms this SOP is followed during internal inspections; initiates CAPA for repeat deviations.
Document Control Officer	Maintains the master copy of this SOP; distributes controlled revisions; retires superseded versions per Section 6.

3 Materials and Equipment

Item	ID / Model	Function
Calibrated refrigerator (2–8 °C)	COLD-04	Short-term reagent storage
Chest freezer (–20 °C)	FRZ-11	Working-stock storage
Ultra-low freezer (–80 °C)	ULT-02	Long-term standard storage
Continuous temperature logger	LOG-TEMP-06	15-min interval monitoring, alarms to QA
Barcode label printer	PRN-BC-3	Internal lot-number labeling
Reagent Receiving Log (Form QF-014-A)	–	Chain-of-custody record on receipt

Never place a reagent into ULT-02 without confirming the logger battery indicator reads green. A dead logger battery was the root cause of a 6-hour undetected excursion in deviation DEV-2025-118.

4 Procedure

4.1 Receiving and Inspection

Step 1

Confirm the shipment matches the packing slip: catalog number, lot number, quantity, and expiration date. Reject and quarantine any shipment with a damaged cold-pack or a temperature-indicator strip showing excursion.

Step 2

Record the reagent name, supplier lot number, receipt date, and assigned internal lot number (format NXB-YYYY-####) in Form QF-014-A within 1 hour of receipt.

Step 3

Print and affix an internal barcode label to every individual container. The label must include the internal lot number, receipt date, and assigned storage condition.

4.2 Storage Assignment

Step 4

Consult the Safety Data Sheet (SDS) storage recommendation and assign the reagent to COLD-04, FRZ-11, or ULT-02 accordingly. Light-sensitive reagents are additionally wrapped in amber sleeves before placement.

Step 5

Place the reagent in the shelf position designated for its category on the storage map posted inside each unit door. First-in-first-out (FIFO) rotation applies within each category.

4.3 Ongoing Monitoring

Step 6

Verify the continuous logger reading against the unit's independent thermometer at the start of each shift and initial Form QF-014-B. Investigate any deviation of more than 2 °C from setpoint immediately.

Step 7

Review the weekly logger export every Monday. Flag any excursion event exceeding 30 minutes to the Laboratory Supervisor for disposition of the affected lots.

A confirmed excursion of more than 2 hours at 2–8 °C, or any duration above –65 °C for ULT-02 contents, requires immediate quarantine of the affected lot and notification of the QA Reviewer before any further use.

5 Quality Control and Deviation Handling

5.1 Acceptance Criteria

- ▶ Storage unit temperature remains within the validated range at every twice-daily check, with no unexplained gap greater than 8 hours between logged readings.
- ▶ 100% of received containers carry a legible internal barcode label before first use.

- ▶ Monthly QA audit finds zero reagents past expiration still present in active storage.

5.2 Deviation Response

Any step in Section 4 that cannot be completed as written (e.g., logger malfunction, mislabeled shipment, unit door left ajar) is a deviation and must be recorded per SOP-QC-002 within 24 hours. The table below summarizes disposition by severity.

Severity	Example	Disposition
Minor	Label printed but placement delayed >1h	Correct and note in log; no lot impact.
Major	Excursion 30 min–2h within 2–8 °C range	Supervisor reviews stability data before re-lease.
Critical	Excursion >2h, or any ULT-02 door-open alarm >15 min	Quarantine lot; QA-led CAPA per SOP-QC-002.

6 Revision History and Approval

6.1 Revision History

Rev.	Date	Description of Change	Author
1.0	2022-03-11	Initial release.	R. Okafor
2.0	2023-07-19	Added ULT-02 procedure; introduced Form QF-014-B.	R. Okafor
3.0	2024-11-04	Revised excursion thresholds after DEV-2024-071 CAPA.	L. Fontaine
4.0	2026-06-15	Added critical-severity quarantine trigger; clarified FIFO rotation mapping.	L. Fontaine

6.2 Approval Signatures

Prepared By

Name

Laila Fontaine

Title

QA Manager

Signature**Date**

Approved By	
Name	Marcus Ideh
Title	Director of Quality Assurance
Signature	
Date	

Next scheduled review: **15 June 2028**. Direct questions regarding this SOP to the Document Control Office at qms@nexabiosciences.example.

References

- [1] International Organization for Standardization. ISO/IEC 17025:2017 – General requirements for the competence of testing and calibration laboratories. Technical report, International Organization for Standardization, Geneva, Switzerland, 2017.